



QUANTORI

BioTechX

# Efficient Data Management for CryoEM HPC in AWS Cloud

# Quantori

**is a Full-Service Scientific Informatics, Data Sciences, and Digital Solutions Provider for the Life Science and Healthcare Industries.**

Using our domain knowledge and technical expertise, we develop cutting-edge data science, digital engineering, and technology platforms for biotech, pharmaceutical, and healthcare companies that accelerate drug discovery and improve patient outcomes.

Our innovative approach harnesses the power of data engineering and informatics, machine learning, emerging technologies, cloud, and HPC expertise to advance research and development and ultimately bridge the gap between meaningful data and patient success.

# CryoEM Processing

Since 2017 CryoEM has become an important tool for drug discovery. The majority of large companies in life science joined the competition then. They built or adjusted existing computation centers to run CryoEM.

CryoEM not only requires having access to an expensive microscope to get RAW data but also infrastructure built for processing and storing the data.

Just entering the field requires significant investment, and small-mid-sized companies naturally select the cloud for CryoEM processing.



# CryoEM Requirements




A single CryoEM run can capture thousands of images and can generate anywhere between 1TB to 15TB of raw datasets



Processing requires extremely fast parallel storage and scalable GPU compute close to the data



Processing of one sample takes days, and in some cases months till achieving a reasonable result



1

RAW data comes from CryoEM LAB into buffer storage:

- Local storage close to microscope (large institutions)
- **Amazon S3 bucket**
- USB drive 😊

2

Data is copied over to shared storage connected to the cluster

3

Processing data is being generated during the run

4


Analysis is completed, data goes to archive (S3 deep glacier)

## Standard Approach

- Persistent GPU instances (or cluster) running 24/7
- Persistent storage used for RAW and intermediate computation data connected to GPU instances

## Advanced Approach

- **Scalable GPU cluster**
- Persistent storage used for RAW and intermediate computation data connected to GPU instances



# Issues with Persistent Storage

1

**Expensive (~\$1000/month for 10TB sample) or slow**

2

**Constantly growing**

3

POSIX filesystem – data is not organized and tagged:

- Multiple copies
- “Trash” data

4

Data catalog - is hard to implement due to the nature of the filesystem


# Dynamic Filesystems Approach in AWS

1 One centralized source – Amazon S3

2 Data can be natively tagged in Amazon S3 as it is object storage

3 Parallel filesystems are being created and removed on-demand

4 Fixed low price for Amazon S3. Fast and expensive filesystems are created on-demand only when needed





+

## Pros

- Pay for what you use – **significant price reduction**
- Data is tagged naturally when enters the system
- Data lifecycle can be applied easily
- Data is isolated on the filesystem level – performance increase

-

## Cons

- Development required to adjust to specific data flow

## Challenge

**The client needed an infrastructure with scalable high-speed storage to transform the Cryo-EM datasets into high-resolution protein structures. For this goal Quantori implemented a virtual GPU/CPU hybrid cluster with a multi-tier storage.**

## Solution

- Hybrid of Scientific platform & Cryo-EM accelerator.
- AWS cloud-based.
- Run Cryo-EM pipeline on multiple samples at one time.
- Integrated RELION and CryoSparc with the HPC scheduler.
- Results are available for post-processing via console/Python Notebooks/ Rstudio.
- LustreFSx and AWS S3 for scratch and long-term storage.

## Benefits

- Platform performance benchmarks fit customer expectation and official RELION/CryoSPARC benchmarks.
- Ability to create/remove Lustre FSx partitions for RAW data significantly reduces storage and platform support costs.
- Research time is reduced due to the availability of results immediately after the computation.
- Batch processing gives a way of the pipeline's steps automation – which increases efficiency and reduces costs.

### Projects/Target Add Dataset +

| Project/Target name                              | Drug      | Environment | Date       |
|--|-----------|-------------|------------|
| ▼ ABCB1 <span style="font-size: 0.8em;">1</span> |           |             |            |
| ▲ ABCB2 <span style="font-size: 0.8em;">3</span> |           |             |            |
|  | • Drug_13 | Env_761     | 01.07.2022 |
|  | • Drug_81 | Env_814     | 12.06.2022 |
|  | • Drug_56 | Env_552     | 07.05.2022 |
| ▼ ABCB3 <span style="font-size: 0.8em;">6</span> |           |             |            |

☰ Dataset summary Uploaded raw data Upload ↑

|   |   |
|---|---|
| <p>Project/Target name:<br/><b>ABCB2</b></p> <p>Drug:<br/><b>Drug_81</b></p> <p>Environment:<br/><b>Env_814</b></p> <p>Date:<br/><b>12.06.22</b></p> <p>S3 bucket:<br/><small>ready@biochem</small></p> <p>Note:<br/>Lorem ipsum dolor sit amet, consectetur adipiscing elit. Aenean commodo ligula eget dolor.</p> | <p>Data source:<br/><b>213123</b> <small>🔗</small></p> <p>Raw pixel size:<br/><b>12</b> <small>🔗</small></p> <p>Accelerating voltage:<br/><b>123123</b> <small>🔗</small></p> <p>Total exposure dose:<br/><b>-</b> <small>🔗</small></p> <p>S3 URL:<br/><small>s3://ready@biochem/ABCB2/Drug_12...</small> <small>🔗</small></p> |
|---|---|

Added Lustras raw\_data\_03 Add Project Lustre +

|   |  |
|---|--|
| <p>Lustre name:<br/><b>raw_data_03</b></p> <p>Type:<br/><b>raw</b></p> <p>S3 URL:<br/><small>s3://ready@biochem/ABCB2/Drug_81/Env_814/120622/project/Target_81</small> <small>🔗</small></p> | <p>Size:<br/><b>2400 MB</b></p> <p>Status:<br/><span style="color: green; font-weight: bold;">● Activated</span></p> |
|---|--|

Deactivate Lustre 🗑

# Testimonial

“ *We are very pleased with the results using Quantori's solution. Being able to process multiple CryoEM datasets in parallel makes handling our constant influx of data possible. We are no longer limited computationally, but as it should be in biochemistry.*

*Whenever we run into an issue or have additional features that we would like to be added, the people at Quantori have been great, and I often get a response from Quantori, if not a solution, in about a day. Overall I recommend Quantori for their knowledgeable, professional staff capable of providing solutions cost-effectively and efficiently.*

**Bharat Reddy**

Senior Scientist

Rectify Pharmaceuticals





[quantori.com](https://quantori.com)